

October 13, 2022

Thermo Scientific 1064Defender Handheld Raman Analyzer

Software Version 1.3.0.2668

RELEASE NOTES

We are pleased to announce the release of the 1064Defender handheld Raman analyzer software version 1.3.0.2668, the latest software for 1064Defender handheld Raman chemical identification systems. This software version adds five new translations and resolves bug fixes.

Release date: September 19, 2022

Version: 1.3.0

Build: 2668

Highlights include:

1) Language Translations

This version introduces support for 5 new languages: Portuguese, Japanese, Korean, Arabic and Polish.

2) Bug Fixes

Several bugs were corrected to improve performance when using the User Library.

- Corrected NO match issue in the User Library
- Resolved export completion issues for User Library items
- Resolved identification issues with User Library Items
- Addressed incorrect duplication of File Names upon export from User Library

For technical support please contact: FSI-Product.Support@thermofisher.com